



CERTIFICATE OF ANALYSIS

Item : GUTTA PERCHA POINTS

1. COMPOSITION

a) Gutta Percha Raw Material	20%
b) Zinc Oxide	60%
c) Barium Sulfate	15%
d) Titanium Dioxide	1%
e) Magnesium Carbonate	2.5%
f) Stearic acid	0.5%
g) Pigment	0.3-1%

2. QUALITY

- a) Stiff but flexible to provide ideal working qualities
- b) Excellent radiopacity
- c) Long shelf life and freshness assured
- d) Largest selection of size and packaging

3. APPEARANCE : Pinkish small cones

4. FLASH POINT AND METHOD OF USE

Gutta Percha Points are used as filling materials in root canal treatment

5. SPECIFICATIONS

- a) Gutta Percha Points are manufactured to comply with American Dental Association specification #57-Type1-Class2
- b) Presented as Gutta Percha Points
Conventional shapes : Accessory sizes
Standardization : I.S.O. Specification

6. LIFE OF PRODUCTS : 3 years from the manufactured date

7. HANDLING : This material is intended for use by properly trained dental professionals only

8. STORAGE : Store in a dry place

Suk Song Oh
Suk Song Oh
META BIOMED CO., LTD.

Suk Song Oh/President
Meta Biomed Co., Ltd.